# beurer

# **MG 16**

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# READ THIS MANUAL COMPLETELY AND CAREFULLY BEFORE USING THIS PRODUCT

Keep this manual in a safe location for future reference

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### 1. IMPORTANT SAFETY NOTES

### Signs and symbols

The following signs appear in the Safety Section and in this manual on pages 3, 4, 5, 6.

READ AND UNDERSTAND THIS ENTIRE MANUAL INCLUDING THE SAFETY SECTION AND ALL INSTRUCTIONS AND WARNINGS COMPLETELY AND CAREFULLY BEFORE USING THIS PRODUCT. FOLLOW ALL SAFETY INSTRUCTIONS AND WARNINGS TO AVOID HAZARDOUS SITUATIONS AND TO MAKE CORRECT USE OF THIS PRODUCT.

<b>A</b>	WARNING	WARNING indicates a hazardous situation which, if not avoided, could result in death or serious injury.
A	CAUTION	CAUTION indicates a hazardous situation which, if not avoided, may result in minor or moderate injury.
	NOTICE	NOTICE addresses practices not related to personal injury, such as product and/or property damage.

# A WARNING

- Packaging materials are a deadly hazard for children and can cause suffocation. Remove all packaging materials immediately and keep them away from children at all times.
- Not for use by children. THIS IS NOT A TOY.

- Swallowing batteries and/or battery fluid can be extremely dangerous.
   Keep the battery and mini massager out of the reach of children and disabled persons. Should any person swallow a battery and/or battery fluid, please call 911 immediately.
- Should battery fluid leak and come into contact with your eyes or skin, immediately rinse with plenty of clean water. Please call 911 immediately.
- · Batteries should not be charged or reactivated by any other means.
- Batteries should not be taken apart, thrown in the fire or short circuited.
- Do not step, jump or sit on the device.
- Protect the mini massager against hard knocks, temperature fluctuations and heat sources which are too close (e.g., stoves, heating radiators).
- Do not drop the mini massager or drop any objects on it.
- This product is not a toy. Keep it out of reach of children.
- People with disabilities, such as impairments, activity limitations, and participation restrictions, or who are physically frail, should be assisted by another person when using this device.
- Do not use this product for any purpose other than massaging human beings.

# A CAUTION:

- This unit is intended for home use only. It is not intended for professional
  use in hospitals or other medical facilities. This unit does not support the
  standards required for professional use.
- If a battery has leaked do not touch the battery fluid. Avoid skin contact (e.g., put on protective gloves) and clean the battery compartment with a dry cloth.
- Do not disassemble this unit. You could get injured.

### The massager:

- · should not be used on children
- should not be used if you experience any abnormalities or injury to the body area to be massaged (e.g. burns, open wounds, scars, rashes or other injuries)
- should not be used if you are experiencing any circulation problems or a headache
- should not be used on the head, bones (e.g. joints, spine) or other sensitive areas of the body
- should not be used after operations
- · should not be used during pregnancy
- should not be used after the consumption of medicines or alcohol (reduced cognition!)
- · should not be used when operating machines
- is for external use only

Use of the device does not replace medical consultation and treatment.

- Those with the following conditions should consult their doctor before using this product:
  - Pregnancy
     Malignant Tumor
  - Diabetes
     Conditions requiring rest

  - · pain of unclarified origin
  - if you are not certain whether the device is suitable for you

- A massage session longer than 15 minutes is not recommended
- Before use, make sure that the device and accessories do not show any visible signs of damage. If you are in any doubt, do not use it and contact your dealer or a specified customer service address
- Do not use any attachments that are not recommended or offered accessories by the manufacturer

### NOTICE:

- Protect the mini massager against hard knocks, moisture, dust, chemicals, major temperature fluctuations and heat sources which are too close (e.g., stoves, heating radiators). This may damage the unit. Please see section 6 "Cleaning, maintenance and disposal" and section 7 "Technical specifications" for details.
- Do not use any aggressive solvents, cleaning agents, detergents or any other strong chemicals to clean the device.
- Never immerse and or spill water or any other liquid onto the mini massager, otherwise liquid will enter the mini massager and cause damage.
- Do not wash the mini massager in a washing machine, dishwasher or tumble dryer.
- Do not place or store appliance where it can fall or be dropped into a tub or sink.
- Do not place or drop into water or any other liquid.
- Always use the size and type of battery indicated.
- · Never use rechargeable batteries. This may damage the unit.
- · Replace weak battery before it discharges completely.
- To avoid draining the battery do not place any items on the mini massager when not in use.

- Leaking batteries may damage the device. If you do not intend to use the mini massager for longer periods, remove the battery from the battery compartment before placing the device in storage.
- Batteries can contain toxins that are harmful to the environment. Always dispose of batteries in accordance with applicable local regulations. Do not dispose of batteries with normal house-hold waste.
- Do not mix different types of batteries together (e.g., alkaline with carbonzinc or old batteries with new ones).
- · Remove used batteries promptly.
- Clean the battery contacts and also those of the device prior to battery installation.
- Observe the local regulations for material disposal. Dispose of the device in accordance with local regulations.
- Do not use mobile telephones, microwaves, and other devices which generate strong electrical or electromagnetic fields near the unit. This may result in incorrect operation of the unit.
- Never attempt to repair the instrument or adjust it yourself. This may damage the unit and impair the functions. If you need to have the unit repaired, please contact our customer service. Please see warranty for service contact. Before submitting any complaint, first check the battery and replace it if necessary.
- Changes or modifications to the device will nullify the user warranty.
- Do not place or store appliance where it can fall or be dropped into a tub or sink.
- Do not place or drop into water or any other liquid.
- Turn off appliance when not in use.

### Content

1. Important Safety Notes 2  A WARNING2	5. Operation10 6. Cleaning, maintenance and
▲ CAUTION 3	disposal1
NOTICE5	Cleaning and maintaining 13
2. Package contents7	Disposal1
3. General notes8	7. Technical specifications 14
4. Batteries 8	8. Warranty 14

### SAVE THESE INSTRUCTIONS FOR FUTURE REFERENCE

If the instruction manual is damaged or if you no longer have the instruction manual in your possession, please contact customer service at 1-800-536-0366 or at info@beurer.com

### Dear Customer,

Thank you for choosing this Beurer product. Our name stands for high-quality, thoroughly tested products for applications in the areas of heat, weight, blood pressure, body temperature, pulse, gentle therapy, massage and air. Congratulations on your purchase of the Beurer mini massager.

Please read these instructions completely and carefully before using this product and keep them for later reference. Be sure to make them accessible to other users and observe the information they contain.

# 2. Package contents

- Mini massager
- 3 x AAA batteries
- . Operating Instructions

### 3. General notes

The device is only intended for the purpose described in these operating instructions. The manufacturer cannot be held liable for damage occurring as a result of inappropriate or careless use.

This device is intended only for the massaging of individual parts of the human body.

The massage must be pleasant and relaxing at all times. Stop the massage or change your position or the contact pressure if you find it is painful or unpleasant.

Use of the device does not replace medical consultation and treatment.

You should consult your doctor before using the massage device if:

- you are fitted with a cardiac pacemaker, implant or other aid,
- in case of pains of unclarified origin,
- if you are not certain whether the device is suitable for you.

# 4. Batteries

# A WARNING:

- Swallowing batteries and/or battery fluid can be extremely dangerous.
   Keep the battery and mini massager out of the reach of children and disabled persons. Should any person swallow a battery and/or battery fluid, please call 911 immediately.
- Should battery fluid leak and come into contact with your eyes or skin, immediately rinse with plenty of clean water. Please call 911 immediately.
- Batteries should not be charged or reactivated by any other means.
- · Batteries should not be taken apart, thrown in the fire or short circuited.

# A CAUTION:

 If a battery has leaked do not touch the battery fluid. Avoid skin contact (e.g., put on protective gloves) and clean the battery compartment with a dry cloth.

### NOTICE:

- Always use the size and type of battery indicated.
- Do not mix different types of batteries together (e.g., alkaline with carbonzinc or old batteries with new ones).
- Never use rechargeable batteries. This may damage the unit.
- Replace weak battery before it discharges completely.
- To avoid draining the battery do not place any items on the mini massager when not in use.
- Leaking batteries may damage the device. If you do not intend to use the mini massager for longer periods, remove the battery from the battery compartment before placing the device in storage.
- Batteries can contain toxins that are harmful to health and the environment.
   Always dispose of batteries in accordance with applicable local regulations. Do not dispose of batteries with normal household waste.

If present, remove the battery insulating strip or remove the battery's protective film and insert the battery according to the polarity. If the mini massager fails to operate, remove the battery completely and reinsert it. Replace the battery cover carefully.

The trial battery supplied with the unit may have a shorter life. Remove battery if you do not intend to use the mini massager for longer periods before placing the device in storage.

9

# 5. Operation

# **A** WARNING:

- Packaging materials are a deadly hazard for children and can cause suffocation. Remove all packaging materials immediately and keep them away from children at all times.
- Not for use by children. THIS IS NOT A TOY.
- Do not step, jump or sit on the device.
- Protect the mini massager against hard knocks, temperature fluctuations and heat sources which are too close (e.g., stoves, heating radiators).
- People with disabilities, such as impairments, activity limitations, and participation restrictions, or who are physically frail, should be assisted by another person when using this mini massager.
- Do not use this product for any purpose other than massaging human beings.

Use of the device does not replace medical consultation and treatment.

- those with the following conditions should consult their doctor before using this product:
  - Pregnancy
     Malignant Tumor
  - Diabetes
     Conditions requiring rest

  - · pain of unclarified origin,
  - if you are not certain whether the device is suitable for you,

- A massage session longer than 15 minutes is not recommended.
- Before use, make sure that the device and accessories do not show any visible signs of damage. If you are in any doubt, do not use it and contact your dealer or a specified customer service address,
- Do not use any attachments that are not recommended or offered accessories by the manufacturer.


#### NOTICE:

- Do not use mobile telephones, microwaves, and other devices which generate strong electrical or electromagnetic fields near the unit. This may result in incorrect operation of the unit.
- Before use, make sure that the device and accessories do not show any visible signs of damage. If you are in any doubt, do not use it and contact your dealer or a specified customer service address.
- Do not use any attachments that are not recommended or offered as accessories by the manufacturer.

### Operation

Inserting or replacing batteries:

- Twist off the battery compartment lid counter-clockwise.
- Insert three alkaline or normal AAA batteries according to the marking in the battery compartment. Always ensure the correct battery polarity!
- Place the battery compartment lid on the battery compartment. Always ensure the correct orientation of the lid.

# 6. Cleaning, maintenance and disposal

### NOTICE:

- Protect the mini massager against hard knocks, moisture, dust, chemicals, major temperature fluctuations and heat sources which are too close (e.g., stoves, heating radiators). This may damage the unit. Please see section 7 "Technical specifications" for details.
- Do not use any aggressive solvents, cleaning agents, detergents or any other strong chemicals to clean the device.
- Never immerse and or spill water or any other liquid onto the mini massager, otherwise liquid will enter the mini massager and cause damage.
- Do not wash the mini massager in a washing machine, dishwasher or tumble dryer.
- Observe the local regulations for material disposal. Dispose of the device in accordance with local regulations.
- Never attempt to repair the instrument or adjust it yourself. This may damage the unit and impair the functions. If you need to have the unit repaired, please contact our customer service. Please see warranty for

service contact. Before submitting any complaint, first check the battery and replace it if necessary.

• Changes or modifications to the device will nullify the user warranty.

### Cleaning and maintaining

To keep your mini massager in the best condition and protect the unit from damage follow these directions. The mini massager should occasionally be cleaned with a damp cloth and do not use aggressive cleaning agents. The unit is not waterproof, do not immerse the unit in water and do not spill water over the unit. If your mini massager accidentally gets wet immediately clean off the water. Store the unit in a safe and dry location. Do not expose the unit to direct sunlight, extreme hot or cold temperatures or humidity. Thoroughly dry any moisture off the unit before storing. Store the unit safely out of the reach of children. Do not subject the unit to strong shocks, such as dropping the unit on the floor. Always use the unit following the instructions provided in this manual.

### Disposal

### NOTICE:

- Batteries should be disposed of separately from household waste.
   Always dispose of batteries as per your local regulations.
- Observe the local regulations for material disposal. Dispose of the device in accordance with local regulations. If you have any questions, please contact the local authorities responsible for waste disposal.

# 7. Technical specifications

Battery type	3 x AAA batteries
Temperature operating range	41°F to 95°F (5°C to 35°C)
Storage temperature range	14°F to 140°F (-10°C to 60°C)
Size	3.9" x 3.9" x 3.5" (100 mm x 100 mm x 90 mm)
Net weight	160 g

Subject to technical modifications.

# 8. Warranty

# Limited Lifetime Warranty For Original Purchaser

Your Beurer Mini massager, Model MG 16, excluding the battery, is warranted to be free from defects in materials and workmanship for the life of the product under normal conditions of intended use and service. This warranty extends only to the original retail purchaser and does not extend to retailers or subsequent owners.

We will, at our option, repair or replace the Beurer Mini massager, Model MG 16, without additional charge, for any part or parts covered by these written warranties. No refunds will be given. Repair or replacement is our only responsibility and your only remedy under this written warranty. If replacement parts for defective materials are not available, Beurer reserves the right to make product substitutions in lieu of repair or replacement.

For warranty service contact our customer service department at 1-800-536-0366 or at info@beurer.com to provide a description of the problem. If the problem is deemed to be within the scope of the limited lifetime warranty, you will be asked to mail the product at your costs in its original package with proof of purchase, your name, address and phone number. If the problem is not deemed to be within the scope of the limited lifetime warranty, we will provide a quotation for repair respectively replacement and return shipping fee.

This warranty does not cover damage caused by misuse or abuse; accident; the attachment of unauthorized accessories; alteration to the product; improper installation; misapplication; lack of reasonable care with respect to the product; unauthorized repairs or modifications; improper use of electrical/power supply; old worn battery; normal wear; loss of power; dropped product; malfunction or damage of an operating part as a result of failure to comply with instructions for use or to provide manufacturer's recommended maintenance; transit damage; theft; neglect; vandalism; or environmental conditions; loss of use during the period the product is at a repair facility or otherwise awaiting parts or repair; or any other conditions whatsoever that are beyond the control of Beurer. This warranty is void if the product is ever used in a commercial or business environment. The maximum liability of Beurer under this warranty is limited to the purchase price actually paid by the customer for the product covered by the warranty, as confirmed by proof of purchase, regardless of the amount of any other direct or indirect damage suffered by the customer.

This warranty is effective only if the product is purchased and operated in the country in which the product is purchased. A product that requires modifications or adaptation to enable it to operate in any other country than the country for which it was designed, manufactured, approved and/or authorized, or repair of products damaged by these modifications is not covered under this warranty.

THE WARRANTY PROVIDED HEREIN SHALL BE THE SOLE AND EXCLUSIVE WARRANTY. ANY IMPLIED WARRANTIES, OBLIGATIONS, OR LIABILITES, INCLUDING BUT NOT LIMITED TO THE IMPLIED WARRANTY OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, ARE LIMITED IN DURATION TO THE DURATION OF THIS APPLICABLE WRITTEN WARRANTY. Some states do not allow limitations on how long an implied warranty lasts, so the above limitations may not apply to you.

IN NO EVENT SHALL BEURER BE LIABLE FOR ANY SPECIAL, INCIDENTAL, INDIRECT OR CONSEQUENTIAL DAMAGES FOR BREACH OF THIS OR ANY OTHER WARRANTY, EXPRESS, IMPLIED OR ANY OTHER THEORY OF LIABILITY, WHATSOEVER. Some states do not allow the exclusion or limitation of special, incidental, or consequential damages, so the above limitation may not apply to you.

Beurer does not authorize anyone, including, but not limited to, retailers, the subsequent consumer purchaser of the product from a retailer or remote purchaser, to obligate Beurer in any way beyond the terms set forth herein.

This warranty does not extend to the purchase of opened, used, repaired, repackaged and/or resealed products, including but not limited to sale of such products on Internet auction sites and/or products by surplus or bulk resellers. Any and all warranties or guarantees shall immediately cease and terminate in connection with any products or parts thereof which are repaired, replaced, altered, or modified, without the prior explicitly written consent of Beurer.

This warranty gives you specific legal rights, and you may also have other rights which may vary from state to state.

For more information regarding our product line in the USA, please visit: www.beurer.com Distributed by: Beurer North America LP Hallandale Beach, FL 33009 www.beurer.com

Questions or comments? Call toll free 1-800-536-0366 or contact info@beurer.com Made in China.

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